

# Alconox, Inc.

Critical cleaning detergents for laboratory, healthcare and industrial applications

30 Glenn Street, White Plains, NY 10603 USA

phone (914) 948-4040 fax (914) 948-4088 e-mail [cleaning@alconox.com](mailto:cleaning@alconox.com) website [www.alconox.com](http://www.alconox.com)

Document #612

## LIQUI-NOX® SAMPLE TRACE ANALYSIS UNDILUTED CONCENTRATE

<u>ELEMENT</u>	<u>DETECTED</u>	<u>DETECTION LIMIT</u>
Chlorine	250 PPM	LESS THAN 10 PPM
Fluorine	NONE	LESS THAN 10 PPM
Bromine	10 PPM	LESS THAN 10 PPM
Iodine	NONE	LESS THAN 10 PPM
Boron	NONE	LESS THAN 10 PPM
Lead	NONE	LESS THAN 10 PPM
Zinc	NONE	LESS THAN 10 PPM
Cadmium	NONE	LESS THAN 50 PPM
Arsenic	NONE	LESS THAN 50 PPM
Antimony	NONE	LESS THAN 50 PPM
Bismuth	NONE	LESS THAN 10 PPM
Tin	NONE	LESS THAN 10 PPM
Sulfur	3000 PPM	LESS THAN 100 PPM
Mercury	NONE	LESS THAN 0.5 PPM
Phosphorus	NONE	LESS THAN 10 PPM
Cobalt	NONE	LESS THAN 10 PPM
Nickel	NONE	LESS THAN 10 PPM
Vanadium	NONE	LESS THAN 10 PPM
Sodium	20,000 ppm	

Note: LIQUI-NOX detergent is recommended for use at 1% dilutions.

