

Alconox, Inc.

Critical Cleaning Experts

30 Glenn St. White Plains NY 10603 USA Tel 914.948.4040 • Fax 914-948-4088 www.alconox.com • cleaning@alconox.com

Report of Inhibitory Residue Test for LIQUINOX®

Study dates: Mar 3- Mar 13, 2020

Analysis Completion Date: March 15, 2020

Sample Identification: M20-0716.11

Master Lot No: L1H9 (master for lots of the form "L, 1-9, A-M, 0-9, and optional A-Z")

Master Lot changes with any significant raw material or manufacturing change. See Certificate of Analysis

for Master Lot reference for any specific sub lot.

TEST METHOD:

Procedure defined and adhered to as described in section 9, 20th Ed Standard Methods for the Examination of Water and Waste Water. Test performed by an independent laboratory.

Detergent: 1% LIQUINOXTM aqueous solution.

- Group A: Washed with LIQUINOXTM, then rinsed 1X DI water
- Group B: Washed with LIQUINOXTM, then rinsed 12X DI water
- Group C: Washed with LIQUINOXTM, no rinsing
- Group D: Sterile Petri dishes

The glassware was sterilized and pour plates were prepared using a known population of *E. aerogenes*

RESULTS:

	Plate #	CFU/mL	CFU/mL	CFU/mL	CFU/mL
		Group A	Group B	Group C	Group D
	Average	67	70	69	62
	<15%				
Ι	Difference	YES	YES	YES	YES

CONCLUSION: As shown in the table, the average CFU/mL differs less than 15% between Groups A, B, C and D. The results indicate that LIQUINOX® detergent does not exhibit any toxic or inhibitory effects.