

Alconox, Inc.

Critical Cleaning Experts
30 Glenn St. White Plains NY 10603 USA
Tel 914.948.4040 • Fax 914-948-4088
www.alconox.com • cleaning@alconox.com

DETONOX®SAMPLE TRACE ANALYSIS UNDILUTED CONCENTRATE

ELEMENT	<u>DETECTED</u>	DETECTION LIMIT	· =
CHLORINE	None	LESS THAN 20	PPM
FLUORINE	None	LESS THAN 20	PPM
BROMINE	None	LESS THAN 20	PPM
IODINE	None	LESS THAN 20	PPM
CHROMIUM	None	LESS THAN 1	PPM
LEAD	None	LESS THAN 0.5	PPM
ANTIMONY	None	LESS THAN 5	PPM
COPPER	None	LESS THAN 1	PPM
ZINC	None	LESS THAN 10	PPM
BISMUTH	None	LESS THAN 1	PPM
TIN	None	LESS THAN 15.0	PPM
CADMIUM	None	LESS THAN 0.5	PPM
ALUMINUM	None	LESS THAN 3.0	PPM
SODIUM	1-3% (w/w)	INGREDIENT	
POTASSIUM	4-7% (w/w)	INGREDIENT	
PHOSPHORUS	2.4% (w/w)	INGREDIENT	
SULFUR	1-3% (w/w)	INGREDIENT	

RECOMMENDED DILUTION IS 1% (1.25 OZ/GL or 10 mL/L **DETONOX** IN WATER). WHEN DILUTED, THE AMOUNT OF THE ABOVE CHEMICAL TRACES WILL OF COURSE BE DILUTED.